


EMnIM - Correlating Electron Microscopy (EM) and Ion Mobility (IM) spectra

Oxford researchers have developed EMnIM, a piece of software capable of accurately calculating Collisional Cross Sections from Electron Microscopy density maps correlated with Ion Mobility data.



To find out more, contact enquiries@innovation.ox.ac.uk and quote the project number.

The above information is provided 'as is' without conditions or warranties. Oxford University Innovation makes no representation and gives no warranty that it is the owner of the intellectual property rights in the technology described. © Oxford University Innovation Ltd 2026, Buxton Court, 3 West Way, Oxford, OX2 0JB.